

## **Carotid Duplex - Image Protocol**

1. Transverse RIGHT CCA Proximal
2. Transverse RIGHT CCA Mid
3. Transverse RIGHT CCA Distal
4. Transverse RIGHT CCA Bulb
5. Transverse RIGHT CCA Bulb with Color
6. Transverse RIGHT CCA Bifurcation
7. Transverse RIGHT CCA Bifurcation with Color
8. Long RIGHT CCA Proximal
9. Long RIGHT CCA Proximal with Color
10. Long RIGHT CCA Proximal with PW Doppler (measure PSV/EDV)
11. Long RIGHT CCA Mid
12. Long RIGHT CCA Mid with Color
13. Long RIGHT CCA Mid with PW Doppler (measure PSV/EDV)
14. Long RIGHT CCA Distal
15. Long RIGHT CCA Distal with Color
16. Long RIGHT CCA Distal with PW Doppler (measure PSV/EDV)
17. Long RIGHT CCA Bulb
18. Long RIGHT CCA Bulb with Color
19. Long RIGHT ICA Proximal
20. Long RIGHT ICA Proximal with Color
21. Long RIGHT ICA Proximal with PW Doppler (measure PSV/EDV)
22. Long RIGHT ICA Mid
23. Long RIGHT ICA Mid with Color
24. Long RIGHT ICA Mid with PW Doppler (measure PSV/EDV)
25. Long RIGHT ICA Distal
26. Long RIGHT ICA Distal with Color
27. Long RIGHT ICA Distal with PW Doppler (measure PSV/EDV)
28. Long RIGHT ECA Proximal
29. Long RIGHT ECA Proximal with Color
30. Long RIGHT ECA Proximal with PW Doppler (measure PSV/EDV)
31. Long RIGHT SUBCLAV with Color
32. Long RIGHT SUBCLAV with PW Doppler (measure PSV/EDV)
33. Long RIGHT Vertebral
34. Long RIGHT Vertebral with and without Color
35. Long RIGHT Vertebral with PW Doppler (measure PSV/EDV)
36. Repeat 1-28 for Left Side